

POSITIVE CONTROL LYSATE FOR ANTI- ARG1 ANTIBODY

SKU: EBL020-ARG1-P1



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 0.3µg/ml staining of Primary antibody in lysates EBL020-ARG1-P1, approx.38kDa band observed. (Calculated MW of 34.7kDa according to Human NP_000036.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

