


POSITIVE CONTROL LYSATE FOR ANTI- BMP1 ANTIBODY

SKU: EBL010-BMP1-S2C2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in lysates EBL010-BMP1-S2C2, Approx.100kDa band observed (calculated MW of 111kDa according to Human NP_006120.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

