

# POSITIVE CONTROL LYSATE FOR ANTI-CAPG ANTIBODY

**SKU:** EBL002-CAPG-P1

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

---

## SPECIFICATIONS

<b>Formulation</b>	5ug/ul protein in RIPA buffer
<b>Unit Size</b>	50 µl
<b>Storage Instructions</b>	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
<b>Usage Summary</b>	Western Blot: 0.3µg/ml staining of Primary antibody in lysate EBL002-CAPG-P1, calculated MW of 38.5kDa according to NP_001738.2 and 36.9kDa according to NP_001243069.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
<b>Shelf Life</b>	3 months
<b>Shipping Instructions</b>	Refrigerated
<b>Application Type</b>	WB

## GALLERY IMAGES

