

POSITIVE CONTROL LYSATE FOR ANTI- CD68 ANTIBODY

SKU: EBL014-CD68-P2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 2µg/ml staining of Primary antibody in lysates EBL014-CD68-P2, Approx.110kDa band observed (calculated MW of 37.4kDa according to NP_001242.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

