

POSITIVE CONTROL LYSATE FOR ANTI- CIRBP ANTIBODY

SKU: EBL007-CIRBP-P1

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 0.03µg/ml staining of Primary antibody in Nuclear lysates EBL007-CIRBP-P1, Approx.19kDa band observed (calculated MW of 18.6kDa according to NP_001271.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

