

POSITIVE CONTROL LYSATE FOR ANTI- DC-SIGN / CD209 ANTIBODY

SKU: EBL008-CD209-P1

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 2µg/ml staining of Primary antibody in lysates EBL008-CD209-P1, Approx.50kDa band observed (calculated MW of 45.8 kDa according to NP_066978.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

