


# POSITIVE CONTROL LYSATE FOR ANTI- GAPDH ANTIBODY

**SKU:** EBL003-GAPDH-S3G2



---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                                                           |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 5ug/ul protein in RIPA buffer                                                                                                                                                                                                                                             |
| <b>Unit Size</b>             | 50 µl                                                                                                                                                                                                                                                                     |
| <b>Storage Instructions</b>  | Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.                                                                                                                                                                                                 |
| <b>Usage Summary</b>         | Western Blot: 0.1µg/ml staining of Primary antibody in lysates EBL003-GAPDH-S3G2, approx. 37kDa band observed.<br>(Calculated MW of 36.1kDa according to Human NP_002037.2).<br>Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K |
| <b>Shelf Life</b>            | 3 months                                                                                                                                                                                                                                                                  |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                                                              |
| <b>Application Type</b>      | WB                                                                                                                                                                                                                                                                        |

## GALLERY IMAGES

