

POSITIVE CONTROL LYSATE FOR ANTI-GAPDH ANTIBODY

SKU: EBL005-GAPDH-S3

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 0.5µg/ml staining of Primary antibody in lysates EBL005-GAPDH-S3, approx. 37kDa band observed. (Calculated MW of 36.1kDa according to Human NP_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

