

POSITIVE CONTROL LYSATE FOR ANTI- GPR94 / TRA1 ANTIBODY

SKU: EBL011-GPR94-G2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 0.03ug/ml staining of Primary antibody in lysates EBL011-GPR94-G2, Approx. 100kDa band observed (calculated MW of 92.5kDa according to Human NP_003290.1 and Mouse NP_035761.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

