

POSITIVE CONTROL LYSATE FOR ANTI- GPX4 ANTIBODY

SKU: EBL008-GPX4-C2

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 2µg/ml staining of Primary antibody in lysates EBL008-GPX4-C2, Approx.19kDa band observed (calculated MW of 22.2kDa according to Human NP_002076.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

