

POSITIVE CONTROL LYSATE FOR ANTI- IMP1/IGF2BP1 ANTIBODY

SKU: EBL019-IGF2BP1-P1



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in nuclear lysates EBL019-IGF2BP1-P1, approx.70kDa band observed. (Calculated MW of 63.5kDa according to Human NP_006537.3 and Mouse NP_034081.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

