

POSITIVE CONTROL LYSATE FOR ANTI- ITGAL ANTIBODY

SKU: EBL024-ITGAL-P1



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 3µg/ml staining of Primary antibody in lysate EBL024-ITGAL-P1, approx.200kDa band observed. (calculated MW of 129kDa according to NP_002200.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

