


POSITIVE CONTROL LYSATE FOR ANTI- MKRN1 ANTIBODY

SKU: EBL022-MKRN1-P1



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 0.5µg/ml staining of Primary antibody in lysates EBL022-MKRN1-P1, approx.55kDa band observed. (Calculated MW of 53.3kDa according to Human NP_038474.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

