

POSITIVE CONTROL LYSATE FOR ANTI- MS2 / ADAM8 / CD156 ANTIBODY

SKU: EBL004-ADAM8-S2G2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in lysates EBL004-ADAM8-S2G2, Approx. 90kDa band observed (calculated MW of 88.8kDa according to NP_001100.3 and 78.8kDa according to NP_001157962.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

