

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

POSITIVE CONTROL LYSATE FOR ANTI- P2RX4 ANTIBODY

SKU: EBL023-P2RX4-S2G2



Telephone: (650) 697-3600



250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa







SPECIFICATIONS

Formulation 5ug/ul protein in RIPA buffer

Unit Size 50 μl

Storage Instructions

Store lysate at -20°C upon arrival. Avoid repeated freeze

/thawed cycles.

Western Blot: $1\mu g/ml$ staining of Primary antibody in lysates

EBL023-P2RX4-S2G2, approx.48kDa band observed. (Calculated

Usage Summary MW of 43.5kDa according to Rat NP_113782.2). Primary

incubation 1 hour at room temperature. Secondary antibody is

EB2ND-002-A at 1:40K.

Shelf Life 3 months **Shipping Instructions** Refrigerated

Application Type WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: <u>(650)</u> 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa