

POSITIVE CONTROL LYSATE FOR ANTI- P2RX4 ANTIBODY

SKU: EBL023-P2RX4-S2G2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in lysates EBL023-P2RX4-S2G2, approx.48kDa band observed. (Calculated MW of 43.5kDa according to Rat NP_113782.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

