

POSITIVE CONTROL LYSATE FOR ANTI- TOM1L1 ANTIBODY

SKU: EBL009-TOM1L1-P1

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in lysates EBL009, approx. 37+50+55kDa bands observed. (Calculated MW of 38.7kDa according to NP_001308102.1, 44.3kDa according to NP_001308103.1 and 53.0kDa according to NP_005477.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

