

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## POSITIVE CONTROL LYSATE FOR ANTI- TOM1L1 ANTIBODY

SKU: EBL009-TOM1L1-P1



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

**Formulation** 5ug/ul protein in RIPA buffer

**Unit Size** 50 μl

**Storage Instructions** Store lysate at -20°C upon arrival. Avoid repeated freeze

/thawed cycles.

Western Blot:  $1\mu g/ml$  staining of Primary antibody in lysates EBL009, approx. 37+50+55kDa bands observed. (Calculated

Usage Summary

MW of 38.7kDa according to NP\_001308102.1, 44.3kDa

according to NP\_001308103.1 and 53.0kDa according to

NP\_005477.2). Primary incubation 1 hour at room temperature.

Secondary antibody is EB2ND-002-A at 1:40K

**Shelf Life** 3 months

**Shipping Instructions** Refrigerated

**Application Type** WB

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa