


RABBIT ANTI-BTAG ANTIBODY

SKU: EB00100



SPECIFICATIONS

Unit Size	100 µg
Storage Instructions	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Usage Summary	The suggested dilutions: ELISA: 1:5000
Immunogen	This antibody is generated from KLH-conjugated synthetic peptide corresponding to 3BTag epitopes (3x QYPALT) derived from major core protein VP7 in bluetongue virus.
Peptide Sequence	CQYPALTQYPALTQYPALT
Purification Method	This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS. Purified polyclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.
Shipping Instructions	Refrigerated
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Application Type	ELISA, WB

GALLERY IMAGES

