

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## RABBIT ANTI-GOAT IGG ANTIBODY AP CONJUGATED

SKU: EB2ND-002-AP

×

## **SPECIFICATIONS**

Formulation Stabiliser: 10mg/ml Bovine Serum Albumin (BSA), IgG and protease-free. Preservative: 0.1% (w/v) Sodium

Azide.

Unit Size 1mg

Storage Instructions

Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of

enzymatic activity. Expiration date is one (1) year from date of opening product.

Usage Summary Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. ELISA:

nary 1:2,000 -1:10,000 WB: 1:500 - 1: 2,500 IHC: 1:200 - 1:1,000

This product was prepared from monospecific rabbit antiserum by immunoaffinity chromatography using goat

Purification Method IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-peroxidase, anti-rabbit serum, goat IgG and goat serum. It was conjugated to Alkaline Phosphatse according to a modified procedure from

Avarameas and Ternyrock, Immunochemistry 32; 1175, 1971.

Shipping Instructions

Refrigerated

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png