

RABBIT ANTI-GOAT IGG ANTIBODY AP CONJUGATED

SKU: EB2ND-002-AP



SPECIFICATIONS

Formulation	Stabiliser: 10mg/ml Bovine Serum Albumin (BSA), IgG and protease-free. Preservative: 0.1% (w/v) Sodium Azide.
Unit Size	1mg
Storage Instructions	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity. Expiration date is one (1) year from date of opening product.
Usage Summary	Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. ELISA: 1:2,000 - 1:10,000 WB: 1:500 - 1: 2,500 IHC: 1:200 - 1:1,000
Purification Method	This product was prepared from monospecific rabbit antiserum by immunoaffinity chromatography using goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-peroxidase, anti-rabbit serum, goat IgG and goat serum. It was conjugated to Alkaline Phosphatase according to a modified procedure from Avarameas and Ternyrock, Immunochemistry 32; 1175, 1971.
Shipping Instructions	Refrigerated
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png