


RABBIT ANTI-PROTEIN C TAG ANTIBODY

SKU: EB00107



SPECIFICATIONS

Formulation	The purified polyclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.
Unit Size	100 µg
Storage Instructions	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Synonym / Alias Names	Protein C; PROC; APC; PC; PROC1; anticoagulant protein C; autoprothrombin IIA; blood coagulation factor XIV; vitamin K-dependent protein C
Usage Summary	The suggested dilutions: ELISA: 1:5000
Immunogen	This antibody is generated from KLH-conjugated synthetic peptide corresponding to the Protein-C epitope CEDQVDPRLIDGK.
Peptide Sequence	CEDQVDPRLIDGK
Purification Method	This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Application Type	ELISA, WB

GALLERY IMAGES

