


RABBIT ANTI-V5 TAG ANTIBODY

SKU: EB00108



SPECIFICATIONS

Unit Size	100 µg
Storage	
Instructions	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Usage	
Summary	The suggested dilutions: ELISA: 1:1000 Western blotting: 1:100
Immunogen	This antibody is generated from KLH-conjugated synthetic peptide corresponding to the epitope GKPIPPELLGLDST of P and V proteins in SV5.
Peptide	
Sequence	GKPIPPELLGLDST
Purification	
Method	This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS. Purified polyclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.
Shipping	
Instructions	Refrigerated
Product	
Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Application	
Type	ELISA, WB

GALLERY IMAGES

