



EBP-PAIR-001-HRP - EBP-PAIR-001-HRP - Immunizing Peptide

Size: 100µg of dried peptide



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

Use

The peptide was used in the production of EB-PAIR-001-HRP – Anti-SARS-CoV2 Spike protein S1 Matched ELISA pair Antibody. It is currently untested in blocking applications.

Peptide Sequence

PROPRIETARY'-C

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.