



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EBP-PAIR-001-HRP - EBP-PAIR-001-HRP - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB-PAIR-001-HRP – Anti-SARS-CoV2 Spike protein S1 Matched ELISA pair Antibody. It is currently untested in blocking applications.

Peptide Sequence

PROPRIETARY'-C

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.