



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EBP05419-M - EBP05419-M - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB05419-M – Goat Anti-AIF1/IBA1 isoform 1 and 3 (Mouse reactivity only) Antibody. It is currently untested in blocking applications.

Peptide Sequence

C-TGPPAKKAISELP,

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.