

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EBP06333 - EBP06333 - Immunizing Peptide

Size: 100µg of dried peptide



#### Use

The peptide was used in the production of EB06333 – Goat Anti-AKAP10 / D-AKAP2 Antibody. It is currently untested in blocking applications.

## **Peptide Sequence**

C-AQYDQPLEKSTKL

## **Supplied As**

100µg of Lyophilised peptide, formulation De-ionized (DI) water. Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

### **Storage**

Shipped at ambient temperature, store at -20°C.