



## **EBP06810 - EBP06810 - Immunizing Peptide**

Size: 100µg of dried peptide



### **UK Office**

#### **Everest Biotech Ltd**

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

### **Use**

The peptide was used in the production of EB06810 – Goat Anti-Mammaglobin A / SCGB2A2 Antibody. It is currently untested in blocking applications.

### **Peptide Sequence**

C-NPQVSKTEYKEL

### **Supplied As**

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

### **Storage**

Shipped at ambient temperature, store at -20°C.