

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EBP09307 - EBP09307 - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB09307 – Goat Anti-Klk6 / neurosin (mouse) Antibody. It is currently tested in blocking applications. Please see antibody <u>datasheet</u> for more details.

Peptide Sequence

KEKPGVYTDVCTHIR

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water. Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.