



EBP10727 - EBP10727 - Immunizing Peptide

Size: 100µg of dried peptide



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

Use

The peptide was used in the production of EB10727 – Goat Anti-Gsta3/Gsta5/Gst Yc2 (rat) Antibody. It is currently untested in blocking applications.

Peptide Sequence

C-EMDPGIVDNFLLK

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.