

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EBP11502 - EBP11502 - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB11502 – Goat Anti-ABHD14B (aa188-200) Antibody. It is currently untested in blocking applications.

Peptide Sequence

HPCYLDKPEEWHT

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water. Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.