



EBP11839 - EBP11839 - Immunizing Peptide

Size: 100µg of dried peptide



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

Use

The peptide was used in the production of EB11839 – Goat Anti-NCU-G1 (mouse aa243-256) Antibody. It is currently untested in blocking applications.

Peptide Sequence

PDCPSVNERN SIDD

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.