



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EBP12320-B - EBP12320-B - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB12320-B – Goat Anti-carbonic anhydrase XII (aa188-199), Biotinylated Antibody. It is currently untested in blocking applications.

Peptide Sequence

CSHLQHVKYKGQE.

Supplied As

100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Reconstitute pellet in distilled water where required. For Blocking studies, reconstitute pellet in 200µl distilled water to obtain a 0.5mg/ml peptide solution.

Storage

Shipped at ambient temperature, store at -20°C.