

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08137 - Goat Anti-14-3-3 tau / YWHAQ Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: YWHAQ, 14-3-3 protein tau, tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide, 14-3-3, 1C5, HS1, 14-3-3 protein

T-cell, 14-3-3 protein theta, protein tau, tyrosine 3/tryptophan 5 -monooxygenase

activation protein, theta polypeptide

Official Symbol: YWHAQ

Accession Number(s): NP_006817.1

Human GeneID(s): 10971

Immunogen

Peptide with sequence C-DDRKQTIDNSQ, from the internal region of the protein sequence according to NP_006817.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to NP_006817.1). Recommended concentration: 0.03-0.1μg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog



EB08137 (0.03µg/ml) staining of Human Brain (Hippocampus) lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.