



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

goat-anti-SST-antibody - Somatostatin proprotein (aa58-70)

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names:

Immunogen

Peptide with sequence CSEPNQTENDALEP, from the Internal region of the protein sequence according to .

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Immunofluorescence: Strong expression of the protein seen in the plasma membrane and endoplasmic reticulum of U2OS and U251 cells. Recommended concentration: 10µg/ml.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig