



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05180 - Goat Anti-ACK1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ACK1, activated p21cdc42Hs kinase, ACK

Official Symbol: TNK2

Accession Number(s): NP_005772

Human GeneID(s): [10188](#)

Immunogen

Peptide with sequence C-HLLGSWGPAAHHR, from the C Terminus of the protein sequence according to NP_005772.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 85-95kDa band seen in Mouse Heart extracts.

Recommended for use at 4-6µg/ml. Overnight incubation is recommended with this antibody.

Has to be failed as SC (120kDa) and J. Biol. Chem, May 2001; 276: 18392 - 18398 (150kDa)

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse