



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05396 - Goat Anti-ADAM15 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ADAM15, ADAM metallopeptidase domain 15 (metargidin), a disintegrin and metalloproteinase domain 15 (metargidin), MDC15, metargidin, a disintegrin and metalloproteinase domain 15, metalloprotease RGD disintegrin protein, metalloproteinase-like, disintegrin-like, and cysteine-rich protein 15

Official Symbol: ADAM15

Accession Number(s): NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1

Human GeneID(s): [8751](#)

Non-Human GeneID(s): 11490 (mouse)

Important Comments: This antibody is expected to recognize isoforms 2, 4, 5, 6 (NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1 resp.)

Immunogen

Peptide with sequence C-APPPPTVSSLYL, from the C Terminus of the protein sequence according to NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 90-95kDa band seen in Daudi lysate. Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human