



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06318 - Goat Anti-AKAP6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: AKAP6, ADAP6, PRKA6, mAkap, ADAP100, AKAP100, KIAA0311, A kinase (PRKA) anchor protein 6, A-kinase anchor protein 100, MGC165020, A-kinase anchor protein 6, protein kinase A anchoring protein 6

Official Symbol: AKAP6

Accession Number(s): NP_004265.3

Human GeneID(s): [9472](#)

Non-Human GeneID(s): 238161 (mouse), 64553 (rat)

Immunogen

Peptide with sequence LTMSVTLSPQRSQ-C, from the N Terminus of the protein sequence according to NP_004265.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: This antibody showed a 190kDa band in lysates of Mouse Heart and of Mouse Parotid (salivary) gland as described by Wu et al, Am J Physiol Cell Physiol. 2010 May;298(5):C1151-8. PMID: 20164376.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow