

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06318 - Goat Anti-AKAP6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: AKAP6, ADAP6, PRKA6, mAKAP, ADAP100, AKAP100, KIAA0311, A kinase (PRKA) anchor protein 6, A-kinase anchor protein 100, MGC165020, A-kinase

anchor protein 6, protein kinase A anchoring protein 6

Official Symbol: AKAP6

Accession Number(s): NP_004265.3

Human GeneID(s): 9472

Non-Human GenelD(s): 238161 (mouse), 64553 (rat)

Immunogen

Peptide with sequence LTMSVTLSPLRSQ-C, from the N Terminus of the protein sequence according to NP_004265.3.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: This antibody showed a 190kDa band in lysates of Mouse Heart and of Mouse Parotid (salivary) gland as described by Wu et al, Am J Physiol Cell Physiol. 2010 May;298(5):C1151-8. PMID: 20164376.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow