

### **International Office**

### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB06772 - Goat Anti-AKT2 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** AKT2, v-akt murine thymoma viral oncogene homolog 2, HGNC:392, PKBBETA, PRKBB, RAC-BETA, Murine thymoma viral (v-akt) homolog-2, rac protein

kinase beta

Official Symbol: AKT2

Accession Number(s): NP\_001617.1; NP\_001229956.1; NP\_001317440.1

Human GeneID(s): 208

Important Comments: This antibody is expected to recognise all reported iasoforms

NP\_001617.1; NP\_001229956.1; NP\_001317440.1.

## Immunogen

Peptide with sequence C-RYDSLGLLELDQRT, from the C Terminus of the protein sequence according to NP\_001617.1; NP\_001229956.1; NP\_001317440.1.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

# **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:16000.

**Western blot:** Approx. 58-60kDa band observed in Human Skeletal Muscle lysates and in Rat Testis lysates, and approx.60kDa in Human Testis lysates, and in preliminary testing of THP-1 cell lysate (calculated MW of 55.8kDa according to Human NP\_001617.1 and 55.5kDa according to Rat NP\_058789.1). Recommended concentration: 0.5-1µg/ml.

Primary incubation 1 hour at room temperature


IHC: Paraffin embedded Human Testis. Recommended concentration: 4-6µg/ml.

**Flow Cytometry:** A customer reported weak staining on THP-1 cells at a concentration of 10ug/ml.


## **Species Reactivity**

Tested: Human, Rat


Expected from sequence similarity: Human, Mouse, Rat, Dog




EB06772 (1μg/ml) staining of Human Skeletal Muscle (A) and (0.75ug/ml) Human Testes (B) lysate (35μg protein in RIPA buffer). Detected by chemiluminescence.



EB06772 (1µg/ml) staining of Rat Testes lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.



EB06772 (4 $\mu$ g/ml) staining of paraffin embedded Human Testis. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB06772 Negative Control showing staining of paraffin embedded Human Testis, with no primary antibody.