

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06279 - Goat Anti-ALKBH Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ALKBH, ABH, alkB, hABH, alkB, alkylation repair homolog (E. coli),

alkylation repair, alkB homolog **Official Symbol:** ALKBH

Accession Number(s): CAA63047.1

Human GeneID(s): 8846

Immunogen

Peptide with sequence HIEQVFSPSASGK-C, from the N Terminus of the protein sequence according to CAA63047.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western Blot: Preliminary experiments gave an approx 25-30kDa band in Hela whole cell lysates at 1.5ug/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 34kDa according to NP_006011. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human