



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08273 - Goat Anti-APOBEC4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: APOBEC4, apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative), C1orf169, MGC26594, OTTHUMP00000060680, putative C->U-editing enzyme APOBEC-4

Official Symbol: APOBEC4

Accession Number(s): NP_982279.1

Human GeneID(s): [403314](#)

Immunogen

Peptide with sequence C-QTKRNPNTKAQE, from the internal region of the protein sequence according to NP_982279.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary testing showed a band at approx 40kDa in Human Testes lysate after 0.03µg/ml antibody staining (calculated MW of 41.6kDa according to NP_982279.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human