



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08273 - Goat Anti-APOBEC4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: APOBEC4, apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative), C1orf169, MGC26594, OTTHUMP00000060680, putative C->U-editing enzyme APOBEC-4

Official Symbol: APOBEC4

Accession Number(s): NP_982279.1

Human GeneID(s): [403314](#)

Immunogen

Peptide with sequence C-QTKRNPNTKAQE, from the internal region of the protein sequence according to NP_982279.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary testing showed a band at approx 40kDa in Human Testes lysate after 0.03µg/ml antibody staining (calculated MW of 41.6kDa according to NP_982279.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human