

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08768 - Goat Anti-APOL4 (C Terminus) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: apolipoprotein L4, APOL4, apolipoprotein L, 4, APOL-IV, APOLIV,

OTTHUMP00000028937, apolipoprotein L-IV

Official Symbol: APOL4

Accession Number(s): NP_085146.2; NP_663693.1

Human GeneID(s): 80832

Important Comments: This antibody is expected to recognize both reported isoforms

(NP_085146.2 and NP_663693.1).

Immunogen

Peptide with sequence C-ENLNELTHIHQSLK, from the C Terminus of the protein sequence according to NP_085146.2; NP_663693.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.


Western blot: Approx. 37kDa band observed in Human Placenta lysates (calculated MW of 38.7kDa according to NP_085146.2). Recommended concentration: 0.3-1μg/ml. Primary incubation was 1 hour. Two additional bands of unknown identity were also consistently observed at 50kDa and at 70kDa. The band at 50kDa was successfully blocked by incubation with the immunising peptide. The other band at 70kDa however was not blocked by the immunizing peptide and it is therefore a non-specific signal.

IHC: Paraffin embedded Human Skeletal Muscle. Recommended concentration: 5µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08768 (0.5 μ g/ml) staining of Human Placenta lysate (35 μ g protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Detected by chemiluminescence.



 $EB08768~(5\mu g/ml)~staining~of~paraffin~embedded~Human~Skeletal~Muscle.~Steamed~antigen~retrieval~with~citrate~buffer~pH~6,~AP-staining.$