

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05121 - Goat Anti-APS Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: APS, adaptor protein with pleckstrin homology and src homology 2

domains

Official Symbol: SH2B2

Accession Number(s): XP_018156

Human GeneID(s): 10603

Immunogen

Peptide with sequence C-GRARAVENQYSFY, from the C Terminus of the protein sequence according to XP_018156.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western blot: Approx 80-85kDa and 40kDa bands seen in 293 cell lysate. Recommended for use at 1-3 µg/ml.

The 80-85kDa is bang on but 40kDa band unexplainable with SC antibody.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human