



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05121 - Goat Anti-APS Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: APS, adaptor protein with pleckstrin homology and src homology 2 domains

Official Symbol: SH2B2

Accession Number(s): XP_018156

Human GeneID(s): [10603](#)

Immunogen

Peptide with sequence C-GRARAVENQYSFY, from the C Terminus of the protein sequence according to XP_018156.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western blot: Approx 80-85kDa and 40kDa bands seen in 293 cell lysate. Recommended for use at 1-3 µg/ml.

The 80-85kDa is bang on but 40kDa band unexplainable with SC antibody.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human