



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05121 - Goat Anti-APS Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: APS, adaptor protein with pleckstrin homology and src homology 2 domains

Official Symbol: SH2B2

Accession Number(s): XP_018156

Human GeneID(s): [10603](#)

Immunogen

Peptide with sequence C-GRARAVENQYSFY, from the C Terminus of the protein sequence according to XP_018156.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western blot: Approx 80-85kDa and 40kDa bands seen in 293 cell lysate. Recommended for use at 1-3 µg/ml.

The 80-85kDa is band on but 40kDa band unexplainable with SC antibody.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human