



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07389 - Goat Anti-Arachidonate 5-lipoxygenase Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** arachidonate 5-lipoxygenase, RP11-67C2.3, ALOX5, ALOX5, 5-LO, 5LPG, LOG5, arachidonic acid 5-lipoxygenase, leukotriene A4 synthase

**Official Symbol:** ALOX5

**Accession Number(s):** NP\_000689.1

**Human GeneID(s):** [240](#)

**Non-Human GeneID(s):** 11689 (mouse), 25290 (rat)

### Immunogen

Peptide with sequence C-ERNKKKQLPYYLSPD, from the C Terminus of the protein sequence according to NP\_000689.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Approx 85kDa band observed in Human Bone Marrow and Rat Stomach lysates (calculated MW of 77.7.0kDa according to rat NP\_036954). Recommended concentration: 0.3-1µg/ml.

### Species Reactivity

**Tested:** Human, Rat

**Expected from sequence similarity:** Human, Mouse, Rat, Dog



EB07389 (0.3 $\mu$ g/ml) staining of Rat Stomach lysate (35 $\mu$ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.