

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB07389 - Goat Anti-Arachidonate 5-lipoxygenase Antibody

Size: 100µg specific antibody in 200µl

### **Target Protein**

Principal Names: arachidonate 5-lipoxygenase, RP11-67C2.3, ALOX5, ALOX5, 5-LO,

5LPG, LOG5, arachidonic acid 5-lipoxygenase, leukotriene A4 synthase

Official Symbol: ALOX5

Accession Number(s): NP\_000689.1

Human GeneID(s): 240

Non-Human GenelD(s): 11689 (mouse), 25290 (rat)

#### **Immunogen**

Peptide with sequence C-ERNKKKQLPYYYLSPD, from the C Terminus of the protein sequence according to NP\_000689.1.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Approx 85kDa band observed in Human Bone Marrow and Rat Stomach lysates (calculated MW of 77.7.0kDa according to rat NP\_036954). Recommended concentration: 0.3-1µg/ml.

# **Species Reactivity**

Tested: Human, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog



EB07389 (0.3 $\mu$ g/ml) staining of Rat Stomach lysate (35 $\mu$ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.