



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05550 - Goat Anti-ARF4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ARF4, ADP-ribosylation factor 4

Official Symbol: ARF4

Accession Number(s): NP_001651

Human GeneID(s): [378](#)

Non-Human GeneID(s): 11843 (mouse), 79120 (rat)

Immunogen

Peptide with sequence GLTISSLFSRL-C, from the N Terminus of the protein sequence according to NP_001651.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 20-22kDa band seen in Mouse Brain lysate.

Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse, Rat